**Uricult® Trio****Instructions for use • English****Intended use**

A dipslide culture method for diagnosing urinary tract infections by demonstrating microbes in urine.

**Principles of the procedure**

The Uricult Trio dipslide system is based on three agar media. One side of the plastic slide is covered with green CLED medium and the other with both reddish-brown MacConkey medium and a colourless E. coli medium. The CLED medium is intended for determining the total bacterial count. The MacConkey medium is intended for detecting gram-negative organisms. This medium contains bile salts that prevent the growth of gram-positive organisms other than enterococci which may grow as minute colonies. The E. coli medium is intended especially for detecting gram-negative, β-glucuronidase-producing organisms which grow as colonies in various shades of brown or grey on this medium. *Escherichia coli* is the most common β-glucuronidase-producing organism encountered in urinary tract infections. Bile salts contained in the E. coli medium prevent the growth of gram-positive organisms on it. Certain yeasts do grow on the E. coli medium, however.

**Reagents****Contents**

Uricult Trio	Cat. No. 68197
Dipslides	10
Patient labels	10
Instructions for use	1

**Storage**

Store Uricult Trio at 7...25°C protected from draught, temperature fluctuations and light sources. Avoid storage near heat-generating appliances.

**Do not allow to freeze.** The expiry date is marked on the box.

**Warnings and precautions**

Uricult Trio is for *in vitro* diagnostic use only.

Do not use the product beyond the expiry date marked on the box. Wear protective clothing and disposable gloves while handling samples or tests, and wash hands thoroughly afterwards.

Do not use the Uricult Trio if you detect discoloration or dehydration of the agar, separation of the growth media from the plastic slide or evidence of bacterial or fungal growth.

Because any colonies growing on Uricult Trio are actual or potential pathogens, do not touch the growth.

**Sample collection and preparation**

Ideally, urine for bacterial culture should remain in the bladder for four hours prior to sampling. Urine samples may be obtained by voiding (clean-voided midstream urine), catheterisation or suprapubic aspiration.

The sample should be applied onto the Uricult Trio slide immediately after collection. The slide should then at once be returned into its protective tube and the cap closed tightly.

If the urine sample needs to be stored prior to inoculation, it should be maintained refrigerated at 2...8°C no longer than 24 hours.

Uricult Trio test results may be affected if the patient has received anti-infective treatment. The test should not be performed until 48 hours after the final dose of medication.

**Procedure**

- Unscrew the slide from the tube without touching the agar surfaces.
- Holding Uricult Trio by the cap, dip the slide into freshly voided midstream urine so that the agar surfaces are totally immersed. If the volume of urine is too small for this, the agar surfaces can be wetted by pouring urine on them, followed by tilting to ensure complete wetting.
- Allow excess urine to drain from the slide.
- Blot the last drops on absorbent paper.
- Screw the slide tightly back into the tube.
- Fill in the patient label and attach it to the tube.
- Place the tube upright in an incubator (36±2°C) for 16–24 hours. The tube may also be sent to a laboratory for incubation.
- To obtain a colony count (CFU/ml), remove the slide from the tube and compare the colony density with the model chart provided in the kit.










**Note:**

- Negative cultures and complicated or catheter-associated UTI samples are recommended to always be incubated for an additional 24 hours to ensure that slow-growing bacteria are detected.
- The inoculated slide may be incubated immediately or stored or transported to a laboratory for incubation and interpretation. Storage or transportation should not exceed 48 hours at 7...25°C, after which Uricult Trio should be incubated at 36±2°C for 16–24 hours. If the slide has been stored or transported for up to 48 hours, only the presence of growth and the colony count should be recorded from it; the colour reaction may be atypical.
- The inoculated slide may be incubated at room temperature for 1–3 days, after which positive cultures may be sent to a specialised laboratory for further investigation<sup>5</sup>. Negative cultures may be incubated for additional 24 hours to detect slow-growing bacteria<sup>6</sup>.

**Quality control**

Quality control tests are performed on each lot of Uricult Trio dipslides at the time of manufacture. Should the user wish to perform his own quality control, the following procedure is recommended:

- Prepare a 10<sup>3</sup>–10<sup>6</sup> bacteria/ml suspension of each of the following bacteria in sterile saline:
  - Staphylococcus aureus* ATCC 25923
  - Escherichia coli* ATCC 25922
  - Proteus mirabilis* ATCC 12453

**Model Chart • Tableau de référence • Standardbildkarte • Tablas de referencia • Tabela de Referência • Avlesningsmall • Tolkningsmall • Mallitaulu****Typical formulation • Formules • Typische Formulierung • Fórmula típica • Formula típica • Formulazione tipica • Τυπική Σύσταση • Typické složení • Sestava gojíšť • Agar indhold • Agar sammensætning • Sammansättning • Koostumus**

CLED medium	MacConkey medium	E. coli medium
Peptone 10.0 g/l	Peptone 20.0 g/l	Peptone 12.0 g/l
Meat extract 3.0 g/l	Lactose 10.0 g/l	MgSO <sub>4</sub> 0.1 g/l
Lactose 10.0 g/l	Neutral red 0.075 g/l	MnCl <sub>2</sub> 0.01 g/l
L-Cystine 0.13 g/l	Bile salts 0.8 g/l	Fe(II) citrate 0.4 g/l
Bromthymol blue 0.03 g/l		Bile salt mixture 2.4 g/l
		8-hydroxyquinoline β-D glucuronide

**Literature • Bibliographie • Literatur • Bibliografia • Referências Bibliográficas • Bibliografia • Βιβλιογραφία • Literatura • Referenser • Litteratur • Kirjallisuus**

- McAllister TA, Arnett GC, Barr W, Kay P: Assessment of plain dipslide quantitation of bacteruria. *Nephron* 11: 111–122, 1973.
- Kass EH: Bacteruria and the diagnosis of infections of the urinary tract. *Archives of Internal Medicine* 100: 709–714, 1957.
- Mackay JP, Sandys GH: Laboratory diagnosis of infections of the urinary tract in general practice by means of a dip-inoculum transport medium. *British Medical Journal* 2: 1286–1288, 1965.
- NCCLS Publication M22-A: Quality Assurance Standards for Commercially Prepared Microbiological Culture Media; Approved standard. Volume 10: 14, 1990.
- Ekman et al.: Streptocult® ja Uricult-levyjen käyttö päivystysaikana.

**Explanation of symbols • Explication des symboles • Zeichenerklärung • Explicación de los símbolos • Explcação de símbolos**

- Spiegazione dei simboli • Επεξήγηση των σύμβολων • Vysvětlivky používající symboly Pojsnila simbolov • Symbolforklaring
- Förklaring på symboler • Förklaring av symboler • Symbolselkylset

IVD	REF	LOT			

In vitro diagnostic medical device	Catalogue number	Batch code	Temperature limitation	Use by	Manufacturer
Dispositif médical de diagnostic <i>in vitro</i>	Référence	Numéro de lot	Limites de température	a utiliser jusqu'à	Fabricant
In-vitro-Diagnostikum	Bestell Nr.	Lot Nr.	Temperaturbegrenzung	Verwendbar bis	Hersteller
Producto sanitario para diagnóstico <i>in vitro</i>	Número de catálogo (nº)	Número de lote	Limitación de temperatura	Fecha de caducidad	Fabricante
Dispositivo medico-diagnóstico <i>in vitro</i>	Código	Número de lote	Límite de temperatura	Utilizar até	Produtor
Ιατρογενοκόπιο βοηθόπιο του χρηματοποίηση για διάγνωση <i>in vitro</i>	Αριθμός καταλόγου	Αριθμός παρτίδων	Πιστοποιοί θερμοκρασίας	Χρησιμοποιήστε έως	Κορυκευτέσται από
Diagnóstico y salud pública instrumental <i>in vitro</i>	Kataložna številka	Číslo šárže	Teplotní omezení	Použítelel do	Výrobce
Diagnóstică și sănătatea publică instrumentă <i>in vitro</i>	Bestellnummer	Serijska stevilka	Teplotné rozmezie	Uporabno do	Proizvajalec
Medicinsk ustyr til in vitro-diagnostikk	Artiklenummer	Lotnummer	Temperaturbegrensning	Utilizzare entro/Scadenza	Fabrikant
Medicinteknisk produkt avsedd för in vitro-diagnostik	Tuotenumero	Batch nr.	Temperaturbegrensning	Brukes innen	Produsent
Medicinteknisk produkt avsedd för in vitro-diagnostik		Eräkoodi	Lämpötilaraja	Används före	Tillverkare

In vitro -diagnostikkalaatko tarkoitettu lääkinnällinen laite	Σ	Protect from draught and temperature fluctuations	Use by	Manufacturer
		Conserver à l'abri des courants d'air et des fluctuations de température	a utiliser jusqu'à	
		Vor Zug und Temperaturschwankungen geschützt lagern	Verwendbar bis	
		Proteger de las corrientes de aire y cambios de temperatura	Fecha de caducidad	
		Proteger das correntes de ar e variações de temperatura	Utilizar até	
		Proteggere da correnti d'aria e variazioni di temperatura	Utilizzare entro/Scadenza	
		Proteger de correntes d'aire y variaciones de temperatura	Použítelel do	
		Proteger de correntes d'aire y variaciones de temperatura	Uporabno do	
		Proteger de correntes d'aire y variaciones de temperatura	Bruges innen	
		Proteger de correntes d'aire y variaciones de temperatura	Används före	
		Proteger de correntes d'aire y variaciones de temperatura	Käytettävä ennen	

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